GENESY 600

ELECTROTHERAPY 4 CHANNELS | 149 PROGRAMS

PROFESSIONAL FOUR-CHANNEL PRODUCT FOR DEMANDING USERS

53 SPORT | 29 FITNESS | 20 PAIN | 12 MICROCURRENTS | 12 SERIAL SEQUENTIAL STIMULATION | 11 REHAB
7 ACTION NOW | 3 DENERVATED | 1 INCONTINENCE | 1 IONOPHORESIS


TECHNICAL FEATURES

Display
Visible area size 2,6”
Channels
4 indipendent (8 electrodes)
Frequency
0,3-150 Hz *
Pulse amplitude
40-450 µs *
Power
0-120 mA per channel
Power supply
Rechargeable batteries
Size
100x160x35 mm
Weight
450 gr

* According to the type of current

EQUIPMENT

- 1 Bag
- 1 GENESY 600 stimulator
- 4 Cables for electrode connection
- 2 Cables for microcurrents and ionophoresis
- 4 Self-adhesive square electrodes
- 4 Self-adhesive rectangular electrodes
- 1 Charger
- 1 Operating manual



MEDICAL DEVICES
CE
0476



- TENS — TENS
- MENS — EMS
- IONO — IONO
- 3S — 3S
- ACTION NOW — ACTION NOW
- RUN TIME — RUN TIME
- LAST 10 — LAST 10
- FAVORITES — FAVORITES
- MULTI USER — MULTI USER
- PRO PROGRAMMABLE — PROGRAMMABLE
- STIM LOCK — STIM LOCK
- 2+2 — FUNCTION 2+2
- MY TRAINER — MY TRAINER
- AUTO STIM — AUTO STIM
- SyS Syncro Stim — SYNCRO STIM
- Work Time — WORK TIME
- RECHARGEABLE BATTERIES

GLOBUS
ITALIAN EXCELLENCE

GENESY 600

ELECTROTHERAPY 4 CHANNELS | 149 PROGRAMS

Rev. En_05_2022

PROGRAMS LIST

REHAB

- Swollen ankles
- Atrophy recovery
- Hemiplegy-upper limbs
- Hemiplegy-lower limbs
- Recovery after ACL surgery
- Functional recovery
- Ankle re-education
- Leg re-education
- Muscle Spasms
- Shoulder subluxation prevention
- Vastus medialis reinforcement
- Motor point pen
- Agonist/Antagonist
- Muscle reinforcement

PAIN

- Conventional antalgic tens
- Modulated antalgic tens
- Low frequency antalgic tens
- Endorphinic tens
- Knee pain
- Menstrual pain
- Post-surgical pain
- Chronic pain
- Shoulder pain (s.h. syndrome)
- Carpal tunnel
- Trapezius pain
- Rotator cuff tendinitis
- Muscle pain
- Sciatica
- Chronic lumbago
- Cervical pain
- Bursitis-tendinitis
- Bone fractures
- Epicondylitis
- Osteoarthritis

IONOPHORESIS

- Ionophoresis

INCONTINENCE

- Mixed incontinence

MICROCURRENTS

- Epicondylitis
- Scapulohumeral periarthritis
- Contusion
- Sciatica
- Lumbago
- Acute pain
- Articular pain
- Stiff neck
- Whiplash
- Shoulder sprain
- Knee sprain
- Tendon inflammation

DENERVATED

- Triangular (low denerv. level)
- Trapezoidal (intermediate denerv. level)
- Rectangular (high denerv. level)

SPORT

- Capillarization
- Warm-up
- Pre-competition warm-up
- Active recovery
- Maximum strength
- Endurance strength
- Explosive strength
- Aerobic endurance
- Reactivity
- Post-competition recovery
- Decontracting
- Hypertrophy

FITNESS

- Firming
- Bio-Pulse firming
- Sculpting
- Bio-Pulse sculpting
- Toning
- Mass building
- Body sculpting
- Definition
- Jogging
- Anaerobic fitness
- Aerobic fitness
- Cramp prevention

SERIAL SEQUENTIAL STIMULATION

The "3S" program list includes 12 parameter combinations.
The "3S" programs are characterized by a delay activation of the channels 3 and 4 compared with the channels 1 and 2..
The Serial Sequential Stimulation permits to stimulate the musculature in kinetic chain thanks to the differentiated activation times of the muscular groups involved.

ACTION NOW

The Action Now program list includes 7 parameter combinations.
The Action Now programs are particularly useful to link and synchronize the electric stimulation with a voluntary action.

DOMINO S.R.L.
Via Vittorio Veneto, 52
31013 Codognè (TV) - ITALY
(+39) 0438 7933
globuscorporation.com

GLOBUS
ITALIAN EXCELLENCE